

NextGen Protocols

TITLE: Immunofluorescent Staining for γ H2AX in Cells

AUTHORS: Marcus Parrish¹, James Evans², Le Ngo¹, and Bevin P. Engelward³

AUTHORS AFFILIATIONS AND TITLES:

¹Graduate Student, Department of Biological Engineering, Massachusetts Institute of Technology, Cambridge, MA


²Research Scientist, Computational and Systems Biology Initiative, Massachusetts Institute of Technology, Cambridge, MA;

³Professor, Department of Biological Engineering, Massachusetts Institute of Technology, Cambridge, MA

MOTIVATION:

- DNA damage is associated with an increased risk of cancer, aging and disease
- γ H2AX is a widely accepted marker for the response to DNA double strand breaks
- This protocol demonstrates a method to stain and quantify the double strand breaks in an *in vitro* system.

EXAMPLE DATA:



γ H2AX Positive MDCK Cell irradiated with 50Gy of Gamma Radiation

PROTOCOL:

Step 1: Gelatinizing the Coverslips

Overview

This process prepares the coverslips to be adherent to the cells. All of this work should be performed in BSL2 Biosafety cabinet

Active Time: 5 minutes

Total Time: 1 hour 35 minutes

Materials:

Falcon 12-well Flat Bottom Culture Plate with lid (Corning 353043)

18 mm #1 Circular Coverslips (Electron Microscopy Sciences 72229-01)

Gelatin from Porcine skin (Sigma G1890)

NextGen Protocols

- 1.1 Mix 5 g of Gelatin with 500 mL of distilled water
- 1.2 Autoclave the 1% Gelatin solution until sterile
- 1.3 Pipet 1 mL of 1% gelatin solution to well in 12-well plate
- 1.4 Place individual coverslip into each well with gelatin (**Video #1**)
- 1.5 Use pipet to press coverslip to bottom of well (**Video #1**)
- 1.6 Rotate plate to ensure gelatin covers entire coverslip (**Video #1**)
- 1.7 Place 12-well plate into 37°C incubator for 1 hour
Tip: Plate can be in incubator for 1-6 hours without effect

Step 2: Plating the Cells

Overview

This step is for seeding the cells prior to beginning the experiment. All of the cell culture in this step should be performed in a BSL2 biosafety cabinet

Active Time: 20 minutes

Total Time: 20 minutes

Materials:

Trypsinized cell solution
Cell specific media

- 2.1 Quantify the number of cells in the cell solution
- 2.2 Centrifuge vial of trypsinized cells at 1000 RPM for 5 minutes at 4°C
- 2.3 Aspirate trypsin solution from the cells
Tip: Be careful not to aspirate the pellet
- 2.4 Resuspend the cells to 1×10^6 Cells/mL
- 2.5 Remove 12-well plate from incubator and aspirate gelatin
- 2.6 Plate 1 mL of cell-specific media to each gelatinized well in 12-well plate
- 2.7 Pipet 0.3×10^6 cells to each well
- 2.8 Rotate plate 360° to ensure complete coverage of area with cells
- 2.9 Place plate in 37°C incubator overnight

Step 3: Experimentation and Fixation

Overview

This process is for performing the actual experimentation and fixation.

Active Time: 20 minutes

Total Time: 20 minutes

Materials:

4% Paraformaldehyde aqueous solution (Electron Microscope Sciences #157-4)
Sterile Dulbecco's Phosphate Buffered Saline (Life Technologies 14190-144)

- 3.1 Perform desired experiment on cells
Note: This experiment can be any procedure that does not require the removal of the cells from the plate.

NextGen Protocols

- 3.2 Upon completion of the experiment, aspirate liquid in wells of 12-well plate
- 3.3 Pipet 0.5 mL of 4% paraformaldehyde to each well
Note: Paraformaldehyde is extremely toxic and should be used inside a ventilated fume hood.
Note: You cannot fix wells independently of one another, because 4% paraformaldehyde on one side of the slide fixes cells on the other side of the plate
- 3.4 Incubate wells at room temperature for 10 min
- 3.5 Removed paraformaldehyde and wash wells twice with PBS
Note: At this point, the cells can be incubated with 1mL PBS and placed in 4°C for two weeks or until further use

Step 4: Immunofluorescent Staining

Overview

This process is for performing the actual permeabilisation, staining and mounting the slides

Active Time: 40 minutes

Total Time: 3 hours + 1 overnight

Materials:

Fine forceps

Triton x-100 (Sigma-Aldrich X100)

Tris Buffered Saline (BioRad 170-6435)

Bovine Serum Albumin (Sigma-Aldrich A3059)

Anti-phospho-Histone H2A.X(Ser139) Antibody, clone JBW301 (Millipore 05-636)

Goat anti-mouse IgG-Fitc (Santa Cruz sc-2010)

26G x ½ Regular Bevel Needles (Becton Dickinson 305111)

1 mL Tuberculin Slip Tip Syringe (Becton Dickinson 309569)

75X25 mm Single End Frosted Slides (VWR 16004-368)

Prolong Gold Antifade with DAPI (Life Technologies P36935)

- 4.1 Dilute 20 µL of Triton X-100 in 10 mL of TBS to produce permeabilization buffer (0.2% Triton X-100 in TBS)
- 4.2 Dissolve 0.1g of BSA in 10 mL TBS to produce blocking buffer (1% BSA in TBS)
Note: This is enough blocking buffer to performs stains on 8 wells. Scale volume accordingly.
- 4.3 Aspirate PBS from wells of the 12 well plate
- 4.4 Incubate each well with 0.5mL of Permeabilization Buffer for 10 min at room temperature
- 4.5 Aspirate Permeabilization Buffer
- 4.6 Incubate each well with 0.5mL of Blocking Buffer for 1 hour at room temperature
- 4.7 Dilute 6 µL of anti-γH2AX Antibody in 3 mL of Blocking Buffer (1:500 Dilution)
- 4.8 Aspirate Blocking Buffer
- 4.9 Pipet 500 µL of Primary Antibody solution to each well in plate
- 4.10 Incubate each plate overnight at 4°C with gentle rocking
- 4.11 Store blocking buffer at 4°C

Next Day:

- 4.12 Aspirate Primary Cell Solutions
- 4.13 Wash all wells with 0.5mL of TBS to each well twice aspirating after each wash
- 4.14 Dilute 6µL Goat- anti-mouse-Fitc in blocking buffer (1:500 Dilution)

NextGen Protocols

- 4.15 Wash all wells with 0.5 mL of TBS twice aspirating after each wash
- 4.16 Bend Syringe needle to a 90 angle on counter with tip facing down (**Video #2**)
- 4.17 Pipet 5 uL of ProLong Antifade onto glass slide
- 4.18 Using Bent Syringe and forceps, pickup coverslip and flip onto slide (**Video #3**)
- 4.19 Image slide as soon as possible
- 4.20 Store slides at -20C until imaged

Step 5: Imaging and Quantification

Overview

This process is to learn how to image and quantify yH2AX pictures

Active Time: 3 hours

Total Time: 3 hours

Materials:

Fluorescent microscope

- 5.1 Place slide on microscope stage and visualize at 10x under DAPI filter
- 5.2 Switch to 60x magnification and focus
- 5.3 Capture images of the slide first under DAPI filter and subsequently under Fitc filter
- 5.4 Using your image software program, create a Dapi/Fitc merged image
- 5.5 Take 5-10 images per slide throughout the entire slide area
Note: To be properly unbiased, search under DAPI filter
- 5.6 Upon capturing images for every slide, have a colleague blind you to the name of every image
- 5.7 Count the number of Dapi positive cells in that image
- 5.8 Count the number of cells with >5 yH2AX foci in that image
- 5.9 Repeat for remaining images.
- 5.10 Unblind yourself to the images names.
- 5.11 Sum the total number of DAPI positive cells and cells with >5 yH2AX foci for each original slide
- 5.12 Finally quantify the percentage yH2AX positive using the equation below:

$$\frac{\text{Total \# of cells with } > 5 \text{ yH2AX foci}}{\text{Total \# of DAPI Positive Cells}} * 100\% \\ = \text{Percetange yH2AX Positive Cells Per Condition}$$